


Adelante[®]

ADVANCED PEEL AWAY INTRODUCER SYSTEM

LATEX-FREE




The Adelante® Peel Away Introducer System is one of the most advanced introducers available for placement of permanent pacing leads and catheters. It offers exceptional insertion, smooth lead movement, high kink resistance and exceptional peeling performance even in complex conditions.

Adelante® Peel Away Introducer System Features:

- Clear numeric print of French sizes
- Ergonomically designed handles
- Color coding of dilator lock and hub to prevent mismatches
- Enhanced hub break force
- Consistent peel through twin channel sheath design
- Smooth dilator/sheath transition
- Smooth guidewire/dilator transition
- Coated sheath for superior lead passage
- Dilator tip indicator
- Compatible guidewire 0.035"/0.89 mm
- FDA clearance and CE certified

Model # Kit	Model # Sheath / Dilator Only	French Size	Sheath Length	Dilator Length
STANDARD				
AP06016B5	AP06016P	6 F	16 cm	19 cm
AP07016B5	AP07016P	7 F	16 cm	19 cm
AP07516B5	AP07516P	7.5 F	16 cm	19 cm
AP08016B5	AP08016P	8 F	16 cm	19 cm
AP08516B5	AP08516P	8.5 F	16 cm	19 cm
AP09016B5	AP09016P	9 F	16 cm	19 cm
AP09516B5	AP09516P	9.5 F	16 cm	19 cm
AP10016B5	AP10016P	10 F	16 cm	19 cm
AP10516B5	AP10516P	10.5 F	16 cm	19 cm
AP11016B5	AP11016P	11 F	16 cm	19 cm
AP11516B5	AP11516P	11.5 F	16 cm	19 cm
AP12016B5	AP12016P	12 F	16 cm	19 cm
AP12516B5	AP12516P	12.5 F	16 cm	19 cm
LONG				
AP06023B5	AP06023P	6 F Long	23 cm	26 cm
AP07023B5	AP07023P	7 F Long	23 cm	26 cm
AP08023B5	AP08023P	8 F Long	23 cm	26 cm
AP09023B5	AP09023P	9 F Long	23 cm	26 cm
AP10023B5	AP10023P	10 F Long	23 cm	26 cm
AP11023B5	AP11023P	11 F Long	23 cm	26 cm
AP12023B5	AP12023P	12 F Long	23 cm	26 cm
EXTRA LONG (Sheath and Dilator only)				
N/A	AP06030P	6 F Extra Long	30 cm	34 cm
	AP07030P	7 F Extra Long	30 cm	34 cm
	AP08030P	8 F Extra Long	30 cm	34 cm
	AP09030P	9 F Extra Long	30 cm	34 cm
	AP10030P	10 F Extra Long	30 cm	34 cm
	AP11030P	11 F Extra Long	30 cm	34 cm
	AP12030P	12 F Extra Long	30 cm	34 cm
	AP06045P	6 F Extra Long	45 cm	49 cm
	AP07045P	7 F Extra Long	45 cm	49 cm
	AP08045P	8 F Extra Long	45 cm	49 cm
	AP09045P	9 F Extra Long	45 cm	49 cm
	AP10045P	10 F Extra Long	45 cm	49 cm
	AP11045P	11 F Extra Long	45 cm	49 cm
	AP12045P	12 F Extra Long	45 cm	49 cm



Adelante® Sterile Kit Contains:

- One Adelante® peel-away sheath
- One dilator
- One guidewire with a flexible straight and J tip with straightener
- One needle with protector
- One syringe

Available in other kit configurations

US patent D576,276. US patent 7,909,798.

For additional information e-mail us at Adelante@oscor.com

Introducers supplied sterile. For single use only. Available bulk upon request.

Standard lengths are noted above. Other lengths available upon request.