


Adelante® *Radial*
INTRODUCER WITH HEMOSTATIC VALVE

LATEX-FREE

The **Adelante® Radial** is the next generation in transradial access sheath introducers, featuring the latest hemostatic valve technology, for minimal blood loss during diagnostic angiography or percutaneous coronary intervention. Hydrophilic kink-resistant coated sheath for easier insertion and removal of device. Provides excellent elongated tapered dilator tip for optimal vessel access.

Adelante® Radial Introducer Set with Hemostatic Valve Features:

- Ergonomically designed hub
- Reliable hemostatic seal to prevent blood loss and air embolism
- Easy French size identification with color coded hub and 3-way stopcock
- Kink resistant and flexible hydrophilic coated sheath for smooth insertion
- Atraumatic tip of dilator provides smooth transition between dilator and sheath during insertion
- Snap-fit dilator with french size indicator
- Dilator with curve memory effect for easy sheath placement
- Excellent visibility under fluoroscopy
- Complete hemostasis during insertion and withdrawal of devices
- Compatible guidewire 0.018"/0.46 mm



Model Number	French Size	Sheath Length	Dilator Length
RADIAL KIT			
ASR04010B5	4 F	10 cm	15 cm
ASR05010B5	5 F	10 cm	15 cm
ASR06010B5	6 F	10 cm	15 cm
ASR07010B5	7 F	10 cm	15 cm
RADIAL SHEATH AND DILATOR SET			
ASR04010P	4 F	10 cm	15 cm
ASR05010P	5 F	10 cm	15 cm
ASR06010P	6 F	10 cm	15 cm
ASR07010P	7 F	10 cm	15 cm

Supplied sterile. For single use only. Available bulk upon request.

Standard lengths are noted above. Other lengths available upon request.


Adelante® Radial Kit (5 per box)

- One Adelante® Radial sheath
- One Adelante® Radial dilator
- One 0.018" guidewire with a flexible tip
- One 21 gauge needle with protector
- One 10 cc syringe

Adelante® Radial Set (10 per box)

- One Adelante® Radial sheath
- One Adelante® Radial dilator

Available in other kit configurations



For additional information, e-mail us at AdelanteRadial@oscor.com