

Adelante[®] SIGMA

INTRODUCER SET WITH HEMOSTATIC VALVE

LATEX-FREE

The **Adelante® SIGMA** is the next generation in vascular access sheath introducers, featuring the latest hemostatic valve technology, for small to large devices.


Adelante® SIGMA Introducer Set with Hemostatic Valve Features:

- Ergonomically designed hub
- Reliable hemostatic seal to prevent blood loss and air embolism
- Easy French size identification with color coded hub and 3-way stopcock
- Kink resistant and flexible hydrophilic coated sheath for smooth insertion
- Atraumatic tip of dilator provides smooth transition between dilator and sheath during insertion
- Adelante **Sigma Sol** version available with a radiopaque atraumatic tip
- Snap-fit dilator with french size indicator
- Dilator with curve memory effect for easy sheath placement
- Excellent visibility under fluoroscopy
- Complete hemostasis during insertion and withdrawal of devices
- Compatible guidewire 0.035"/0.89 mm


Model Number Standard Tip	Model Number Sol Tip	French Size	Sheath Length	Dilator Length
STANDARD				
ASG04010	ASGS04010	4 F	10 cm	15 cm
ASG05010	ASGS05010	5 F	10 cm	15 cm
ASG06010	ASGS06010	6 F	10 cm	15 cm
ASG07010	ASGS07010	7 F	10 cm	15 cm
ASG08010	ASGS08010	8 F	10 cm	15 cm
ASG09010	ASGS09010	9 F	10 cm	15 cm
ASG10010	ASGS10010	10 F	10 cm	15 cm
ASG11010	ASGS11010	11 F	10 cm	15 cm
ASG12010	ASGS12010	12 F	10 cm	15 cm
LONG				
ASG04023	ASGS04023	4 F Long	23 cm	28 cm
ASG05023	ASGS05023	5 F Long	23 cm	28 cm
ASG06023	ASGS06023	6 F Long	23 cm	28 cm
ASG07023	ASGS07023	7 F Long	23 cm	28 cm
ASG08023	ASGS08023	8 F Long	23 cm	28 cm
ASG09023	ASGS09023	9 F Long	23 cm	28 cm
ASG10023	ASGS10023	10 F Long	23 cm	28 cm
ASG11023	ASGS11023	11 F Long	23 cm	28 cm
ASG12023	ASGS12023	12 F Long	23 cm	28 cm

Introducer kits supplied sterile. For single use only. Available bulk upon request.

Standard lengths are noted above. Other lengths available upon request.



For additional information, e-mail us at AdelanteSigma@oscor.com



Adelante® SIGMA Sterile Kit Contains:

- One Adelante® Sigma sheath
- One dilator
- One guidewire with a flexible straight, J tip and depth markings
- One needle with protector
- One syringe

Available in other kit configurations