


ANGIOS[®] classic

DIAGNOSTIC ANGIOGRAPHY CATHETERS

Features:



Angios® Classic Diagnostic Angiography Catheters:

This angiography catheter features a proprietary HD braid reinforced shaft offering maximum pushability while retaining excellent curve stability, maximum flow rate during contrast injection, and an ultra soft radiopaque catheter tip to minimize vessel trauma. Oscor's Angios® diagnostic catheters are available in a variety of curve, size and length configurations. Other custom convenience kit configurations available upon request.

Angios® Classic Features:

- **Kink resistant** shaft featuring superior **pushability**, optimum **torque** and longitudinal **stability**
- **Hydrophobic Coated** shaft for easy insertion
- **Large** shaft lumen for **high flow rate**
- **Soft atraumatic tip** to minimize vessel trauma
- **Curve retention** for accurate catheter placement
- **Excellent** shaft and tip **visualization** under fluoroscopy
- Ergonomically designed and **color coded hub**
- Ergonomically designed tip straightener for easy loading over guidewire and sheath introducer (Pigtail catheter model only)
- All models available in a wide variety of lengths and curve sizes
- Accepts maximum guidewire up to 0.038"/0.96 mm

Angios® Classic Specifications:

Angios® catheters are available in the following packaging configurations:

Oscor's Angiography Catheter - 100 Series, consists of 100 single catheter pouches per box.

Oscor's Angiography Catheter Kit - 300 Series, consists of 50 single catheter kit pouches per box. Angiography Catheter Kit contains a Judkins left, a Judkins right, and either a Straight Pigtail or a Pigtail 145° in a variety of length configurations.


Oscor's Angiography Convenience Kit - 400 Series, consists of 1 single convenience kit per box. Angiography Convenience Kit contains a 3 catheter kit, a fixed or movable core guidewire, and an Adelante® Sigma introducer kit.

Oscor's Angiography Catheter - 500 Series, consists of 5 single catheter pouches per box.

Custom packaging and curve configurations are available upon request.

For a complete list of models e-mail us at sales@oscor.com

Angios® Classic Inner-lumen Sizes and Lengths:



Products not available in all regions.



P: 727.937.2511
F: 727. 934.9835
E: sales@oscor.com
www.oscor.com

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

