


Curve	Length	French Size	100 Series Model # (100 per box)	500 Series Model # (5 per box)
AR MOD I	100 cm	● 5 F	ANG105-0009	ANG505-0009
		● 6 F	ANG106-0011	ANG506-0011
	110 cm	● 5 F	ANG105-0010	ANG505-0010
		● 6 F	ANG106-0012	ANG506-0012
AR MOD II	100 cm	● 5 F	ANG105-0021	ANG505-0021
		● 6 F	ANG106-0023	ANG506-0023
	110 cm	● 5 F	ANG105-0022	ANG505-0022
		● 6 F	ANG106-0024	ANG506-0024
AL I	100 cm	● 5 F	ANG105-0103	ANG505-0103
		● 6 F	ANG106-0105	ANG506-0105
	110 cm	● 5 F	ANG105-0104	ANG505-0104
		● 6 F	ANG106-0106	ANG506-0106
AL II	100 cm	● 5 F	ANG105-0109	ANG505-0109
		● 6 F	ANG106-0111	ANG506-0111
	110 cm	● 5 F	ANG105-0110	ANG505-0110
		● 6 F	ANG106-0112	ANG506-0112
IM	100 cm	● 5 F	ANG105-0203	ANG505-0203
		● 6 F	ANG106-0205	ANG506-0205
	110 cm	● 5 F	ANG105-0204	ANG505-0204
		● 6 F	ANG106-0206	ANG506-0206
MPA I	100 cm	● 5 F	ANG105-0303	ANG505-0303
		● 6 F	ANG106-0305	ANG506-0305
	110 cm	● 5 F	ANG105-0304	ANG505-0304
		● 6 F	ANG106-0306	ANG506-0306
NIH	100 cm	● 5 F	ANG105-0403	ANG505-0403
		● 6 F	ANG106-0405	ANG506-0405
	110 cm	● 5 F	ANG105-0404	ANG505-0404
		● 6 F	ANG106-0406	ANG506-0406
TALON 1	100 cm	● 5 F	ANG1-0958	ANG5-0958
		● 6 F	ANG1-0961	ANG5-0961
	110 cm	● 5 F	ANG1-0959	ANG5-0959
		● 6 F	ANG1-0962	ANG5-0962
	120 cm	● 5 F	ANG1-0960	ANG5-0960
		● 6 F	ANG1-0963	ANG5-0963

100/500 SERIES SPECIALTY CURVES



5 F 0.069" / 1.75 mm	6 F 0.082" / 2.08 mm
1200 PSI PSI	100 - 120 cm

100 SERIES: 100 single catheter pouches provided per box. 500 SERIES: 5 single catheter pouches provided per box.

Maximum Recommended Guidewire 0.038". Catheter supplied sterile. For single use only.

Standard lengths are noted above. Other lengths available upon request.

