




VEIN PICKS

The Vein Pick is a disposable plastic accessory designed to facilitate introduction of lead, catheter or devices alike; in the event of cut-down procedure.



VEIN PICK – SHORT

Model #	Color	Overall Length	Height at Tip	QTY
VP52Y	● Yellow	45 mm	9 mm	50/ case



VEIN PICK – LONG

Model #	Color	Overall Length	Height at Tip	QTY
VP51B	● Blue	110 mm	10 mm	50/ case



Vein picks are supplied sterile. For single use only. Available bulk upon request.

Products not available in all regions.



P: 727.937.2511
 F: 727. 934.9835
 E: sales@oscor.com
 www.oscor.com

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

