



Adelante®



# Adelante®

ADVANCED PEEL AWAY INTRODUCER SYSTEM

## Features:




### Adelante® Kit Contains:

- One Adelante® peel-away sheath
- One dilator
- One guidewire (up to 0.038"/0.96 mm) with a flexible straight and J tip with straightener
- One needle with protector
- One syringe

### Adelante® Set Contains:

- One Adelante® peel-away sheath
- One dilator

Available in other kit configurations



The Adelante® Peel Away Introducer System is one of the most advanced introducers available for placement of permanent pacing leads and catheters. It offers exceptional insertion, smooth lead movement, high kink resistance and exceptional peeling performance even in complex conditions.

### Adelante® Peel Away Introducer System Features:

- Clear numeric print of French sizes
- Ergonomically designed handles
- Color coding of dilator lock and hub to prevent mismatches
- Optimized hub break force
- Consistent peel through twin channel sheath design
- Smooth dilator/sheath transition
- Smooth guidewire/dilator transition
- Coated sheath for superior lead passage
- Dilator tip indicator
- Accepts maximum guidewire up to 0.038"/0.96 mm

Model # Kit (Five per box)	Model # Set (Ten per box)	French Size	Sheath Length	Dilator Length
<b>STANDARD</b>				
AP06016B5	AP06016P	6 F	16 cm	19 cm
AP07016B5	AP07016P	7 F	16 cm	19 cm
AP07516B5	AP07516P	7.5 F	16 cm	19 cm
AP08016B5	AP08016P	8 F	16 cm	19 cm
AP08516B5	AP08516P	8.5 F	16 cm	19 cm
AP09016B5	AP09016P	9 F	16 cm	19 cm
AP09516B5	AP09516P	9.5 F	16 cm	19 cm
AP10016B5	AP10016P	10 F	16 cm	19 cm
AP10516B5	AP10516P	10.5 F	16 cm	19 cm
AP11016B5	AP11016P	11 F	16 cm	19 cm
AP11516B5	AP11516P	11.5 F	16 cm	19 cm
AP12016B5	AP12016P	12 F	16 cm	19 cm
AP12516B5	AP12516P	12.5 F	16 cm	19 cm
<b>LONG</b>				
AP06023B5	AP06023P	6 F	23 cm	26 cm
AP07023B5	AP07023P	7 F	23 cm	26 cm
AP08023B5	AP08023P	8 F	23 cm	26 cm
AP09023B5	AP09023P	9 F	23 cm	26 cm
AP10023B5	AP10023P	10 F	23 cm	26 cm
AP11023B5	AP11023P	11 F	23 cm	26 cm
AP12023B5	AP12023P	12 F	23 cm	26 cm
<b>EXTRA LONG</b>				
N/A	AP06030P	6 F	30 cm	34 cm
N/A	AP08030P	8 F	30 cm	34 cm
N/A	AP10030P	10 F	30 cm	34 cm
N/A	AP12030P	12 F	30 cm	34 cm
N/A	AP07045P	7 F	45 cm	49 cm
N/A	AP08045P	8 F	45 cm	49 cm
N/A	AP09045P	9 F	45 cm	49 cm

Introducers supplied sterile. For single use only. Available bulk upon request. Standard lengths are noted above. Other lengths available upon request.

For additional information e-mail us at [sales@oscor.com](mailto:sales@oscor.com)

**Products not available in all regions.**



P: 727.937.2511  
F: 727.934.9835  
E: [sales@oscor.com](mailto:sales@oscor.com)  
[www.oscor.com](http://www.oscor.com)

**CAUTION:** Federal Law restricts this device to sale by or on the order of a physician (U.S.A). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

