


# Breezeway II

## DELIVERY SHEATH

Transseptal Guiding Sheath with optional Loader

The Breezeway® II Delivery Sheath is intended for the introduction of diagnostic and therapeutic devices into the human vasculature, including intracardiac.


## SHEATH DIAMETERS



## LOADER DIAMETERS



# ORDERING INFORMATION



### Kit Includes:


- (1) Breezeway II Sheath
- (1) Breezeway II Dilator

### Plus Kit Includes:

- (1) Breezeway II Sheath
- (1) Breezeway II Dilator
- (1) Loader


## BREEZEWAY II – 6F - 8F

Model # Kit <i>(One per box)</i>	Model # Plus Kit <i>(One per box)</i>	A1 Sheath Inner Diameter	A2 Usable Sheath Length	A3 Sheath Tip Curve	B1 Usable Loader Length
BG060079045	BG060079045L	● 6 F	79 cm	45°	8.5 cm
BG060110045	BG060110045L	● 6 F	110 cm	45°	8.5 cm
BG070079045	BG070079045L	● 7 F	79 cm	45°	8.5 cm
BG070079180	BG070079180L	● 7 F	79 cm	180°	8.5 cm
BG070110045	BG070110045L	● 7 F	110 cm	45°	8.5 cm
BG080079045	BG080079045L	● 8 F	79 cm	45°	12.5 cm
BG080079180	BG080079180L	● 8 F	79 cm	180°	12.5 cm
BG080110045	BG080110045L	● 8 F	110 cm	45°	12.5 cm



## BREEZEWAY II – 9F - 14F

Model # Kit <i>(One per box)</i>	Model # Plus Kit <i>(One per box)</i>	A1 Sheath Inner Diameter	A2 Usable Sheath Length	A3 Sheath Tip Curve	B1 Usable Loader Length
BG090079045	BG090079045L	● 9 F	79 cm	45°	12.5 cm
BG090079180	BG090079180L	● 9 F	79 cm	180°	12.5 cm
BG090110045	BG090110045L	● 9 F	110 cm	45°	12.5 cm
BG100079045	BG100079045L	● 10 F	79 cm	45°	12.5 cm
BG100079180	BG100079180L	● 10 F	79 cm	180°	12.5 cm
BG100110045	BG100110045L	● 10 F	110 cm	45°	12.5 cm
BG110079045	BG110079045L	● 11 F	79 cm	45°	16.5 cm
BG110079180	BG110079180L	● 11 F	79 cm	180°	16.5 cm
BG120079045	BG120079045L	● 12 F	79 cm	45°	16.5 cm
BG140079045	BG140079045L	● 14 F	79 cm	45°	16.5 cm



Delivery sheath supplied sterile. For single use only. Available bulk upon request.  
Standard lengths are noted above. Other lengths and curves available upon request.

Products not available in all regions.



P: 727.937.2511  
F: 727.934.9835  
E: sales@oscor.com  
www.oscor.com

**CAUTION:** Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

