



Adelante® VBT

The Adelante® Valve Bypass Tool is designed to be used to temporarily open (dilate) a hemostatic introducer valve for ease of lead passage and protecting the lead integrity. The VBT should be used anytime a defibrillation lead coils or porous tip pass through a hemostatic introducer valve (e.g., at implant or repositioning).

Adelante® Valve Bypass Tool Features:

- Ergonomically designed handles
- Color coding of dilator lock and hub to prevent mismatches
- Enhanced hub break force and easy peeling
- Packaging: Box of 25/single units

ADELANTE VBT — 3 CM SERIES

Model/Catalog #	French Size	Sheath Length
A06VBT03B25	● 6 F	3 cm
A07VBT03B25	● 7 F	3 cm
A75VBT03B25	● 7.5 F	3 cm
A08VBT03B25	● 8 F	3 cm
A85VBT03B25	● 8.5 F	3 cm
A09VBT03B25	● 9 F	3 cm
A95VBT03B25	● 9.5 F	3 cm
A10VBT03B25	● 10 F	3 cm
A105VBT03B25	● 10.5 F	3 cm
A11VBT03B25	● 11 F	3 cm
A12VBT03B25	● 12 F	3 cm
A125VBT03B25	● 12.5 F	3 cm

ADELANTE VBT — 5 CM SERIES

Model/Catalog #	French Size	Sheath Length
A06VBT05B25	● 6 F	5 cm
A07VBT05B25	● 7 F	5 cm
A75VBT05B25	● 7.5 F	5 cm
A08VBT05B25	● 8 F	5 cm
A85VBT05B25	● 8.5 F	5 cm
A09VBT05B25	● 9 F	5 cm
A95VBT05B25	● 9.5 F	5 cm
A10VBT05B25	● 10 F	5 cm
A105VBT05B25	● 10.5 F	5 cm
A11VBT05B25	● 11 F	5 cm
A12VBT05B25	● 12 F	5 cm
A125VBT05B25	● 12.5 F	5 cm

VBT supplied sterile. For single use only. Available bulk upon request.

Standard lengths are noted above. Other lengths available upon request.

Products not available in all regions.

